

PREPARE TO COLLECT YOUR WATER SAMPLES

- 1 Choose one faucet you're going to use for all water sample collection – kitchen, bathroom, outside.
- 2 Wash your hands prior to collection.
- 3 Put on the enclosed nitrile gloves.
- 4 Remove any attachments, aerators and screens on your faucet.
- 5 Sanitize the end of the faucet with the enclosed alcohol swab, or with household bleach on a cloth (no more than 5 seconds).
- 6 Run cold water for 5 minutes at a typical slow rate that you would use to fill a glass of water to drink.
- 7 Do not touch the inside of the cap or container, and do not rinse out any container. Several contain a preservative necessary for the test.